

# Why the LOCalizer™ system may be the right choice for you.



## About Hologic

**Hologic empowers people to live healthier lives everywhere, every day.**

We achieve this by bringing **The Science of Sure®** to life. The Science of Sure reflects who we are as a leading innovator in women's health, helping healthcare professionals around the world diagnose and treat their patients with precision, certainty and confidence.

With groundbreaking technology at the core, our innovations are designed to achieve exceptional clinical results, making it possible to detect, diagnose and treat illnesses and other health conditions earlier and more effectively.

### Convenient scheduling

Your doctor can schedule for the RFID Tag to be implanted into the breast anytime prior to surgery.

### Immediate confirmation

The surgeon can confirm the removal of the precise Tag during surgery by the Tag's unique ID.

### The path towards better outcomes

The LOCalizer system may potentially reduce the chance of follow-up surgery.



Guide to your  
**LOCalizer™ wire-free  
guidance system**  
procedure

To learn more about LOCalizer™ and Hologic, visit us at **Hologic.com**

MISC-05501 Rev.001 (4/19) All rights reserved. Hologic, The Science of Sure and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. Hologic is an exclusive distributor and licensee of the LOCalizer product and trademark. The LOCalizer product is manufactured by Health Beacons. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.

**LOCalizer™**  
Wire-free Guidance System

# A guide to breast surgeries

# How the LOCalizer™ system works

# What to expect with the LOCalizer™ system

**Due to imaging and screening advancements, breast cancer is often caught at earlier stages. However, early detection means that lesions are typically smaller and therefore harder to locate during surgery. Today, about 50% of all lesions are “non-palpable,” which means your surgeon cannot find them by touch alone.**

Traditional breast lesion localization methods involved using surgical wires or radioactive seeds as markers, but these techniques may potentially have disadvantages for doctors and patients. Often involving inconvenience and multiple procedures on the same day, these procedures can possibly lead to unnecessary stress and anxiety.

Now there is a way for your healthcare team to mark and find breast lesions— with the LOCalizer™ that potentially enables precision and ease of use for breast surgery guidance.



Introducing LOCalizer™, the wireless breast lesion localization system from Hologic.

The LOCalizer™ system is a wireless radio frequency identification (RFID) breast lesion localization system that potentially enables precision and ease of use for breast surgery guidance.

The LOCalizer™ system features:

- **Advanced Miniature RFID Tag** designed to be implanted into the breast at any time prior to the day of surgery
- **Real-Time Distance Measurement** reports distance measurement to the Tag.
- **Pencil-sized single-use Surgical Probe** guides the surgeon towards the Tag during the operation.
- **Visible RFID Tag** can be placed under and is visible by ultrasound and x-ray.

## How does the Tag placement work?

The LOCalizer™ system marks the lesion with a miniature radiofrequency Tag designed to be implanted into the breast anytime prior to the day of surgery. Each Tag has a unique identification number and is implanted with a needle, guided with an ultrasound or x-ray machine. Your doctor will use the LOCalizer™ system reader to confirm the position and ID of the Tag. This information may help your surgeon plan the surgical path and removal of the lesion. Once the Tag has been implanted you can return home and resume normal activity.

## What happens during surgery?

Your surgeon will determine the location of the Tag with the LOCalizer™ system reader and probe. The handheld reader will display the distance from the Tag throughout the surgery. By measuring the distance to the target lesion, the LOCalizer™ system may facilitate accurate surgical excision.

Once the lesion has been removed, your surgeon will use LOCalizer™ system probe and reader to confirm that the Tag has been taken out with the tissue.